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
Customer ID 394721

Please scan QR Code on the LFD strip

Product Name RIDA QUICK Aflatoxin RQS

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





# RIDA® SMART APP – on-site mycotoxin testing with smartphone-based test evaluation

Michael Fischer, David Steinmann, Verena Zimmermann, Walter Lübke and Michael Mättner, R-Biopharm AG, Darmstadt, Germany

## Introduction

Mycotoxin contaminations of food and feed have a huge economic impact and pose a serious risk to human and animal health. Mycotoxin contaminations of crops are unavoidable but mycotoxins can be managed. Monitoring mycotoxin contaminations by testing is necessary to verify that the products will meet international regulations and guidelines. Yet, instead of testing large numbers of end-products, a more

pro-active approach would have many benefits. During the entire process from field to food or feed critical steps can be identified to monitor mycotoxins. For this approach a mobile, easy to use and quick tool for on-site mycotoxin testing is the optimal solution. Lateral flow-based tests are well accepted methods for this. However, for a quantitative result, a lab environment is generally required.

## Aim of this study

R-Biopharm AG now presents the next generation in rapid on-site mycotoxin testing, the RIDA® SMART APP. This software application allows you to use a smartphone as a powerful and easy to handle reader for a flexible and fast mycotoxin testing with RIDA® QUICK lateral flow tests – anytime, anywhere. It scans the test strip and

automatically calculates the results by analyzing the line pattern on the taken picture. Here we show data that demonstrates the accuracy and robustness of the RIDA® SMART APP software. R-Biopharm RIDA® QUICK mycotoxin tests can be evaluated reliably and fast with this innovative and highly flexible tool.

## Material and methods

RIDA® QUICK DON A.L. No. R5908	Trilogy® reference material (wheat)							
Target value	ND	0.5	0.9	1.6	2.1	3.5	4.5	6.2
Recovery [%]	-	106	93	93	100	95	104	93
RIDA® QUICK Aflatoxin RQS A.L. No. R5205	Trilogy® reference material (corn)							
Target value	ND	0.5	1.1	1.9	2.7	3.6	4.8	6.2
Recovery [%]	-	113	104	102	105	100	104	97
RIDA® QUICK Aflatoxin RQS ECO A.L. No. R5206	Trilogy® reference material (corn)							
Target value	ND	1.7	5.9	14.3	20.2	31.6	50.8	98.7
Recovery [%]	-	-	97	110	104	106	96	82
RIDA® QUICK Zearalenon RQS A.L. No. R5504	Trilogy® reference material (corn)							
Target value	ND	1.7	5.9	14.3	20.2	31.6	50.8	98.7
Recovery [%]	-	-	88	84	93	88	100	97
RIDA® QUICK Fumonisin RQS A.L. No. R5606	Trilogy® reference material (corn)							
Target value	ND	0.6	1.0	2.2	3.2	6.8	9.2	12.5
Recovery [%]	-	92	111	98	94	101	84	-
RIDA® QUICK T-2 / HT-2 RQS A.L. No. R5304	Oats sample (for spike level see target values)							
Target value	ND	50	100	200	400	600	800	1000
Recovery [%]	-	103	113	96	93	93	92	91

Fig. 1: Recovery of RIDA® QUICK mycotoxin test strips analyzed with RIDA® SMART APP

We analyzed Trilogy® reference material with all available R-Biopharm AG quantitative RIDA® QUICK mycotoxin lateral flow tests and evaluated the results with the RIDA® SMART APP software installed on commercially available NEXUS 6 smartphones. Here we show representative recovery data of all RIDA® QUICK mycotoxin tests as well as information on the minimal influence of

surrounding light conditions on the test results. We also demonstrate the consistency of results of the RIDA® SMART APP software when it is installed on NEXUS 6 smartphones from different sources. For this we purchased nine phones from around the world and compared results that were generated by using a reference test strip.

R-Biopharm quantitative RIDA® QUICK tests were performed according to instructions for use and evaluated with RIDA® SMART APP. As samples naturally contaminated Trilogy® reference materials were used. Certified contamination values were set as target values. For the analysis of oats with RIDA® QUICK T-2 / HT-2 RQS, spiking experiments were performed since no Trilogy® reference materials were available. Non-detect (ND) samples were

defined as such by Trilogy® (by HPLC or LC-MS/MS values below the detection limit). Recovery values are given as (%) if samples tested with a RIDA® QUICK test and evaluated with RIDA® SMART APP gave results below the detection limit or above the upper end of the measurement range. Representative example data is shown in the table below as recovery values, given in percent (%).


Fig. 2: Influence of different light conditions

To validate influence of different light conditions on RIDA® SMART APP results Trilogy® reference materials were tested according to instructions for use and evaluated with RIDA® SMART APP installed on a Nexus 6 smartphone. The graph shows results of aflatoxin ND (non-detect) samples and samples contaminated with Aflatoxin (max. level 75 µg/kg) analyzed with a commercially available lot of R-Biopharm RIDA® QUICK Aflatoxin RQS lateral flow test. Values were generated by using one NEXUS 6 phone. The same extracts were used at different light conditions. Light conditions were as follows: a) light box 19V, b) light box 19V, c) lab bench dark (closed blinds), d) lab bench sunny (open blinds, sunlight), e) lab bench rainy (open blinds, rainy day). The individual lines show the consistency of the 5 recorded curves under different light conditions and controlled as well as non-controlled environments. X-axis shows the target values of the reference materials. Y-axis shows the ratio of the color intensity of test line to control line.


Fig. 3: Inter reader variation

To validate the inter reader variation of NEXUS 6 smartphones nine individual phones (reader) have been tested on the same reference strip (tests strip showing a consistent line pattern) and under the same conditions. Results are given as mg/kg. Standard deviations (SD) and the coefficients of variation (CV) of the measured values were calculated (n=20) for each NEXUS 6. The coefficients of variation were between 3.52 % and 5.43 % and the results were very consistent. All measurements were performed under controlled conditions (Light Box 19V).

## Conclusion

We were not only able to demonstrate that the innovative evaluation tool RIDA® SMART APP gives accurate and precise results under controlled light conditions by showing very good recoveries with Trilogy® reference material. We also show that the RIDA® SMART APP software is a reliable tool to work in a non-controlled environment. The software generates reliable results in combination with all quantitative R-Biopharm RIDA® QUICK mycotoxin tests available on the US market. Additional to that powerful features turn this software into a fast, valuable and cost-effective alternative to conventional lateral flow readers. Results are automatically stored in the internal application database and reports can be recalled any time. By using the connectivity options of the smartphone results can be exported via Wi-Fi or Bluetooth to connected network printers. Optional, results can be sent by E-Mail from any production site worldwide. A Nexus Smartphone with installed RIDA® SMART APP software is a highly mobile lateral flow evaluation device that is independent from any power supply. This makes it an optimal tool for testing in combination with lateral flow mycotoxin tests – mobile, fast, accurate.

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Google Inc. Nexus 6, Nexus 6P